

Benchtop Centrifuge with Swing Buckets, 4000 rpm

Cat. No. Q5120

Buy now for only **\$1500 USD**



The ideal low speed centrifuge for general lab applications

- Compact design saves you valuable lab space
- Sturdy, stainless steel chamber provides reliable performance

Versatile and flexible to suit a variety of experimental needs

- Speed range of 100-4,000 rpm, adjustable in increments of 100 rpm
- Comes with rotor containing: **4 Swing Buckets (8 x 15ml + 4 x 50 ml)**
- Optional accessories available separately:

2 Swing Buckets (4 x 50 ml)

2 Swing Buckets (8 x 15 ml)

2 x 2 x 96 ELISA Plate

Simplify routine operations with easy-to-use controls

- Turn a knob to quickly change parameters
- Push a button to easily convert between RPM and RCF



Goes up to 4,000 rpm
or 2,720 x g



Can hold 15 & 50 ml
tubes or ELISA plate



1 Year Manufacturer's
Warranty

SPECIFICATIONS

Default Configuration 8 x 15 ml + 4 x 50 ml tubes

Optional Configurations: 16 x 15 ml tubes

8 x 50 ml tubes

2 x 2 x 96 ELISA plate

Max Speed/RCF..... 4,000 rpm / 2,720 x g

Max Accel/Decel Time 10 s / 15 s

Time Range 1-99 minutes

Noise Level..... ≤60 dB(A)

Power Supply 120 V

Dimensions (D x W x H) 280 x 500 x 320 mm

Total Weight..... 35 kg

Auto Lid Release Yes, with manual release cord

Safety Alarms Automatic imbalance detection

Certifications..... ISO 13485:2003, CE, REACH, RoHS